

Emerin Rabbit mAb


ABZ5138

Product Overview

Name	Emerin Rabbit mAb
Accession(Primary)	P50402
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 29 kDa; Observed MW: 29 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:500; IF-1:20-1:50; FC-1:20-1:50

Target information (P50402)

Gene Name	EMD
About Gene Symbol	EMD (Emerin) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Emerin in 293T lysates
using Emerin antibody.