


HP1 gamma Rabbit mAb

ABZ5139

Product Overview

Name	HP1 gamma Rabbit mAb
Accession(Primary)	Q13185
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 21 kDa; Observed MW: 24 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:2000

Target information (Q13185)

Gene Name CBX3**About Gene Symbol**

Western blot analysis of HP1 gamma in (1) HeLa lysates; (2) NIH/3T3 lysates using HP1 gamma antibody.

