


Slit2 Rabbit mAb

ABZ5141

Product Overview

Name	Slit2 Rabbit mAb
Accession(Primary)	O94813
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 170 kDa; Observed MW: 58 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (O94813)

Gene Name SLIT2**About Gene Symbol**

Western blot analysis of Slit2 in SH-SY5Y lysates using Slit2 antibody.

