


GABARAPL2/GATE-16 Rabbit mAb

ABZ5142

Product Overview

Name	GABARAPL2/GATE-16 Rabbit mAb
Accession(Primary)	P60520
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P60520)

Gene Name GABARAPL2**About Gene Symbol**

Western blot analysis of GABARAPL2 in (1) HeLa lysates; (2) mouse spleen lysates; (3) rat brain lysates using GABARAPL2 antibody.

