


MOG Rabbit mAb

ABZ5143

Product Overview

Name	MOG Rabbit mAb
Accession(Primary)	Q16653
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 28 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:5000

Target information (Q16653)

Gene Name MOG**About Gene Symbol**

Western blot analysis of Myelin oligodendrocyte glycoprotein in mouse brain lysates using MOG antibody.

