

Insulin Receptor Related Protein Rabbit mAb

ABZ5146

Product Overview

Name	Insulin Receptor Related Protein Rabbit mAb
Accession(Primary)	P14616
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 144 kDa; Observed MW: 80 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (P14616)

Gene Name	INSRR
About Gene Symbol	