


PKC eta Rabbit mAb

ABZ5148

Product Overview

Name	PKC eta Rabbit mAb
Accession(Primary)	P24723
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 78 kDa; Observed MW: 40,78 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (P24723)

Gene Name PRKCH**About Gene Symbol**

Western blot analysis of PKC eta in MCF-7 lysates using PKC eta antibody.

