


OTX1/2 Rabbit mAb

ABZ5149

Product Overview

Name	OTX1/2 Rabbit mAb
Accession(Primary)	P32242
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 37 kDa; Observed MW: 32 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; IP-1:20

Target information (P32242)

Gene Name OTX1/OTX2**About Gene Symbol**

Western blot analysis of Otx1/2 in Y79 lysates
using OTX1/2 antibody.

