


RAP80 Rabbit mAb

ABZ5153

Product Overview

Name	RAP80 Rabbit mAb
Accession(Primary)	Q96RL1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 80 kDa; Observed MW: 80 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-P-1:200-1:500

Target information (Q96RL1)

Gene Name UIMC1**About Gene Symbol**

Western blot analysis of RAP80 in 293T lysates using RAP80 antibody.

