

PCK2 Rabbit mAb


ABZ5158

Product Overview

Name	PCK2 Rabbit mAb
Accession(Primary)	Q16822
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 71 kDa; Observed MW: 71 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:200-1:500; FC-1:20-1:50

Target information (Q16822)

Gene Name	PCK2
About Gene Symbol	



Western blot analysis of PCK2 in (1) HepG2 lysates; (2) mouse brain lysates; (3) rat spleen lysates using PCK2 antibody.