

SF2 Rabbit mAb


ABZ5161

Product Overview

Name	SF2 Rabbit mAb
Accession(Primary)	Q07955
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 28 kDa; Observed MW: 33 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:100-1:500; FC-1:50-1:100

Target information (Q07955)

Gene Name	SRSF1
About Gene Symbol	



Western blot analysis of SF2 in mouse spleen lysates using SF2 antibody.