

mtTFA Rabbit mAb


ABZ5162

Product Overview

Name	mtTFA Rabbit mAb
Accession(Primary)	Q00059
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 29 kDa; Observed MW: 24 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; IP-1:10-1:100

Target information (Q00059)

Gene Name	TFAM
About Gene Symbol	



Western blot analysis of mtTFA in MCF-7 lysates
using mtTFA antibody.