


Glutamate Receptor 2/3 Rabbit mAb

ABZ5169

Product Overview

Name	Glutamate Receptor 2/3 Rabbit mAb
Accession(Primary)	P42262
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 99 kDa; Observed MW: 99 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P42262)

Gene Name GRIA2/GRIA3**About Gene Symbol**

Western blot analysis of GluR2+GluR3 in (1) Human fetal brain lysates; (2) mouse brain lysates using Glutamate Receptor 2/3 antibody.

