

VPS35 Rabbit mAb


ABZ5170

Product Overview

Name	VPS35 Rabbit mAb
Accession(Primary)	Q96QK1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 92 kDa; Observed MW: 92 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:200-1:2000; FC-1:10-1:100

Target information (Q96QK1)

Gene Name	VPS35
About Gene Symbol	VPS35 (VPS35 Retromer Complex Component) is a protein kinase and signaling regulator involved in vesicular trafficking. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of VPS35 in (1) HeLa lysates; (2) mouse kidney lysates using VPS35 antibody.