

ALDH2 Rabbit mAb


ABZ5171

Product Overview

Name	ALDH2 Rabbit mAb
Accession(Primary)	P05091
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 56 kDa; Observed MW: 56 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:2000; IP-1:10-1:100

Target information (P05091)

Gene Name	ALDH2
About Gene Symbol	



Western blot analysis of ALDH2 in HepG2 lysates
using ALDH2 antibody.