

SATB2 Rabbit mAb

ABZ5174


Product Overview

Name	SATB2 Rabbit mAb
Accession(Primary)	Q9UPW6
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 83 kDa; Observed MW: 83 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (Q9UPW6)

Gene Name SATB2

About Gene Symbol



Western blot analysis of SATB2 in HT-1080 lysates using SATB2 antibody.

