

TdT Rabbit mAb

ABZ5176


Product Overview

Name	TdT Rabbit mAb
Accession(Primary)	P04053
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 59 kDa; Observed MW: 59 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (P04053)

Gene Name

DNTT

About Gene Symbol

Western blot analysis of TdT in Jurkat lysates
using DNA Nucleotidylexotransferase antibody.

