

GRK3 Rabbit mAb


ABZ5178

Product Overview

Name	GRK3 Rabbit mAb
Accession(Primary)	P35626
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 80 kDa; Observed MW: 80 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:20-1:50

Target information (P35626)

Gene Name	GRK3
About Gene Symbol	



Western blot analysis of GRK3 in Ramos lysates
using GRK3 antibody.