

alpha Internexin Rabbit mAb


ABZ5179

Product Overview

Name	alpha Internexin Rabbit mAb
Accession(Primary)	Q16352
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 55 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:20-1:50; IF-1:20-1:50; IF-1:20; IP-1:20

Target information (Q16352)

Gene Name	INA
About Gene Symbol	



Western blot analysis of alpha Internexin in SH-SY5Y lysates using alpha Internexin antibody.