

SFPQ Rabbit mAb


ABZ5180

Product Overview

Name	SFPQ Rabbit mAb
Accession(Primary)	P23246
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 76 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:1000-1:5000; ICC/IF-1:50-1:100; FC-1:20-1:50

Target information (P23246)

Gene Name	SFPQ
About Gene Symbol	SFPQ (Splicing Factor Proline And Glutamine Rich) is a protein phosphatase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms



Western blot analysis of SFPQ in A431 lysates
using SFPQ antibody.