

Geminin Rabbit mAb


ABZ5182

Product Overview

Name	Geminin Rabbit mAb
Accession(Primary)	O75496
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 24 kDa; Observed MW: 30 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200

Target information (O75496)

Gene Name	GMNN
About Gene Symbol	GMNN (Geminin DNA Replication Inhibitor) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Geminin in HeLa lysates
using Geminin antibody.