

KDM4A Rabbit mAb


ABZ5184

Product Overview

Name	KDM4A Rabbit mAb
Accession(Primary)	O75164
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 121 kDa; Observed MW: 150 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; IP-1:50-1:100

Target information (O75164)

Gene Name	KDM4A
About Gene Symbol	



Western blot analysis of KDM4A / JHDM3A / JMJD2A in (1) HEK293 lysates; (2) NIH3T3 lysates using KDM4A antibody.