

ZYX Rabbit mAb


ABZ5187

Product Overview

Name	ZYX Rabbit mAb
Accession(Primary)	Q15942
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 61 kDa; Observed MW: 78 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:200-1:500; FC-1:50-1:100; IP-1:20-1:50

Target information (Q15942)

Gene Name	ZYX
About Gene Symbol	



Western blot analysis of Zyxin in MCF-7 lysates
using ZYX antibody.