


Human IgD Rabbit mAb

ABZ5188

Product Overview

Name	Human IgD Rabbit mAb
Accession(Primary)	P01880
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 42 kDa; Observed MW: 84 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:2000-1:5000

Target information (P01880)

Gene Name IGHD**About Gene Symbol**

Western blot analysis of IGHD in Human tonsil lysates using Human IgD antibody.

