


Sodium/Hydrogen Exchanger 1/NHE-1 Rabbit mAb

ABZ5189

Product Overview

Name	Sodium/Hydrogen Exchanger 1/NHE-1 Rabbit mAb
Accession(Primary)	P19634
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 91 kDa; Observed MW: 130 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P19634)

Gene Name SLC9A1**About Gene Symbol**

Western blot analysis of Sodium / Hydrogen Exchanger 1 in MOLT4 lysates using SLC9A1 antibody.

