


Prostate Secretory Protein Rabbit mAb

ABZ5190

Product Overview

Name	Prostate Secretory Protein Rabbit mAb
Accession(Primary)	P08118
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 13 kDa; Observed MW: 13 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000

Target information (P08118)

Gene Name MSMB**About Gene Symbol**

Western blot analysis of Prostate Secretory Protein/PSP in Human prostate cancer lysates using Prostate Secretory Protein antibody.

