


SMARCA5 Rabbit mAb

ABZ5192

Product Overview

Name	SMARCA5 Rabbit mAb
Accession(Primary)	O60264
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 122 kDa; Observed MW: 122 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200

Target information (O60264)

Gene Name SMARCA5**About Gene Symbol**

Western blot analysis of SNF2H in MCF-7 lysates using SMARCA5 antibody.

