


TSH beta Rabbit mAb

ABZ5193

Product Overview

Name	TSH beta Rabbit mAb
Accession(Primary)	P01222
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 16 kDa; Observed MW: 16 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500

Target information (P01222)

Gene Name TSHB**About Gene Symbol**

Western blot analysis of TSH beta in Human pituitary lysates using Thyrotropin Subunit beta antibody.

