

KLF4 Rabbit mAb


ABZ5198

Product Overview

Name	KLF4 Rabbit mAb
Accession(Primary)	O43474
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; ICC/IF-1:50-1:100; FC-1:20-1:50; IP-1:20-1:50

Target information (O43474)

Gene Name	KLF4
About Gene Symbol	KLF4 (KLF Transcription Factor 4) is a transcription factor involved in angiogenesis. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of KLF4 in HeLa lysates
using KLF4 antibody.