


CCR8 Rabbit mAb

ABZ5201

Product Overview

Name	CCR8 Rabbit mAb
Accession(Primary)	P51685
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:100; FC-1:500-1:2000

Target information (P51685)

Gene Name CCR8**About Gene Symbol**

Western blot analysis of CCR8 in rat large intestine lysates using CCR8 antibody.

