


Integrin beta 4 Rabbit mAb

ABZ5202

Product Overview

Name	Integrin beta 4 Rabbit mAb
Accession(Primary)	A2A863
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 202 kDa; Observed MW: 202 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:100-1:200

Target information (A2A863)

Gene Name Itgb4**About Gene Symbol**

Western blot analysis of Integrin beta 4 in SW480 lysates using Integrin beta 4 antibody.

