

Apolipoprotein L1/APOL1 Rabbit mAb


ABZ5203

Product Overview

Name	Apolipoprotein L1/APOL1 Rabbit mAb
Accession(Primary)	O14791
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 44 kDa; Observed MW: 44 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:20-1:50; FC-1:10-1:100; IP-1:10-1:100

Target information (O14791)

Gene Name	APOL1
About Gene Symbol	



Western blot analysis of Apolipoprotein L1/APOL1
in Human serum lysates using ApoL1 antibody.