


KA1 Rabbit mAb

ABZ5205

Product Overview

Name	KA1 Rabbit mAb
Accession(Primary)	Q16099
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 107 kDa; Observed MW: 107 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q16099)

Gene Name GRIK4**About Gene Symbol**

Western blot analysis of KA1 in mouse brain
lysates using KA1 antibody.

