


CD14 Rabbit mAb

ABZ5206

Product Overview

Name	CD14 Rabbit mAb
Accession(Primary)	P08571
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 40 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200

Target information (P08571)

Gene Name CD14**About Gene Symbol**

Western blot analysis of CD14 in Human tonsil
lysates using CD14 antibody.

