

LI Cadherin Rabbit mAb

ABZ5207


Product Overview

Name	LI Cadherin Rabbit mAb
Accession(Primary)	Q12864
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 92 kDa; Observed MW: 120 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500

Target information (Q12864)

Gene Name CDH17

About Gene Symbol



Western blot analysis of LI Cadherin in HCT116, SW480 lysates using LI Cadherin antibody.

