

MDB3 Rabbit mAb


ABZ5208

Product Overview


Name	MDB3 Rabbit mAb
Accession(Primary)	O95983
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 33 kDa; Observed MW: 33 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:1000; FC-1:10-1:100; IP-1:20

Target information (O95983)

Gene Name	MBD3
About Gene Symbol	



Western blot analysis of MDB3 in C6 lysates using MDB3 antibody.



Western blot analysis of MDB3 in A375, A431, mouse brain lysates using MDB3 antibody.