

Tuberin Rabbit mAb


ABZ5209

Product Overview

Name	Tuberin Rabbit mAb
Accession(Primary)	P49815
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 201 kDa; Observed MW: 201 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; FC-1:20-1:50

Target information (P49815)

Gene Name	TSC2
About Gene Symbol	TSC2 (Tuberin) is a protein phosphatase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Tuberin in U87-MG
lysates using Tuberin antibody.