

STAT2 Rabbit mAb


ABZ5210

Product Overview

Name	STAT2 Rabbit mAb
Accession(Primary)	P52630
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 98 kDa; Observed MW: 98 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:20-1:50; FC-1:50-1:100

Target information (P52630)

Gene Name	STAT2
About Gene Symbol	



Western blot analysis of STAT2 in K562 lysates
using STAT2 antibody.