

Filamin A Rabbit mAb


ABZ5211

Product Overview


Name	Filamin A Rabbit mAb
Accession(Primary)	P21333
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 281 kDa; Observed MW: 281 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:50-1:100; FC-1:50-1:100

Target information (P21333)

Gene Name	FLNA
About Gene Symbol	



Western blot analysis of Filamin A in rat brain lysates using Filamin A antibody.



Western blot analysis of Filamin A in HepG2, 3T3, Hela lysates using Filamin A antibody.