

CD146 Rabbit mAb


ABZ5212

Product Overview

Name	CD146 Rabbit mAb
Accession(Primary)	P43121
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 72 kDa; Observed MW: 120 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P43121)

Gene Name	MCAM
About Gene Symbol	MCAM (Melanoma Cell Adhesion Molecule) is a cell surface marker involved in adaptive immunity and angiogenesis. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of CD146 in HeLa, A375, Raji lysates using CD146 antibody.