


PLC gamma 1 Rabbit mAb

ABZ5217

Product Overview

Name	PLC gamma 1 Rabbit mAb
Accession(Primary)	P19174
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 149 kDa; Observed MW: 149 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P19174)

Gene Name PLCG1**About Gene Symbol**

Western blot analysis of PLC gamma 1 in Jurkat lysates using PLCG1 antibody.

