

VCAM1 Rabbit mAb


ABZ5219

Product Overview


Name	VCAM1 Rabbit mAb
Accession(Primary)	P19320
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 81 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:100-1:200; FC-1:20-1:50; IP-1:20-1:50

Target information (P19320)


Gene Name	VCAM1
About Gene Symbol	



Western blot analysis of VCAM1 in mouse kidney lysates using VCAM1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using VCAM1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of VCAM1 in K562 using VCAM1 antibody.

