

Ubiquitin Rabbit mAb


ABZ5220

Product Overview


Name	Ubiquitin Rabbit mAb
Accession(Primary)	P0CG47
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 26 kDa; Observed MW: 8 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:100; FC-1:50-1:100

Target information (P0CG47)


Gene Name	UBB
About Gene Symbol	



Western blot analysis of Ubiquitin in 293T lysates using Ubiquitin antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon using Ubiquitin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of Ubiquitin in Raji using Ubiquitin antibody.

