

Cytokeratin 9 Rabbit mAb


ABZ5221

Product Overview

Name	Cytokeratin 9 Rabbit mAb
Accession(Primary)	P35527
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 62 kDa; Observed MW: 62 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; IP-1:10-1:100

Target information (P35527)

Gene Name	KRT9
About Gene Symbol	



Western blot analysis of Cytokeratin 9 in HeLa
lysates using Cytokeratin 9 antibody.