

Transglutaminase 2 Rabbit mAb


ABZ5224

Product Overview

Name	Transglutaminase 2 Rabbit mAb
Accession(Primary)	P21980
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 77 kDa; Observed MW: 77 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:2000-1:5000

Target information (P21980)

Gene Name	TGM2
About Gene Symbol	TGM2 (Transglutaminase 2) is a protease enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic



Western blot analysis of TGM2 in HUVEC lysates
using Transglutaminase 2 antibody.