

Peroxiredoxin 3/PRDX3 Rabbit mAb


ABZ5226

Product Overview

Name	Peroxiredoxin 3/PRDX3 Rabbit mAb
Accession(Primary)	P30048
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 28 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; FC-1:200-1:1000; IP-1:100-1:500

Target information (P30048)

Gene Name	PRDX3
About Gene Symbol	PRDX3 (Peroxiredoxin 3) is a peroxiredoxin antioxidant enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Peroxiredoxin 3 in K562, Raji, MCF-7 lysates using Peroxiredoxin 3 antibody.