


CD204 Rabbit mAb

ABZ5227

Product Overview

Name	CD204 Rabbit mAb
Accession(Primary)	P21757
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 50 kDa; Observed MW: 75 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P21757)

Gene Name MSR1**About Gene Symbol**

Western blot analysis of MSR1 in (1) HeLa lysates;
(2) mouse placenta lysates using MSR1 antibody.

