

CD79a Rabbit mAb


ABZ5228

Product Overview

Name	CD79a Rabbit mAb
Accession(Primary)	P11912
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 25 kDa; Observed MW: 45-55 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500

Target information (P11912)

Gene Name	CD79A
About Gene Symbol	



Western blot analysis of CD79a in Human spleen lysates using CD79a antibody.