

TNF Receptor 2 Rabbit mAb


ABZ5230

Product Overview


Name	TNF Receptor 2 Rabbit mAb
Accession(Primary)	P20333
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 48 kDa; Observed MW: 65 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:20-1:50; FC-1:20-1:100; IP-1:10-1:100

Target information (P20333)


Gene Name	TNFRSF1B
About Gene Symbol	




Western blot analysis of TNF Receptor II in 3T3, C6 lysates using TNF Receptor 2 antibody.



Western blot analysis of TNF Receptor II in Jurkat, HeLa lysates using TNF Receptor II antibody.



Immunocytochemistry analysis of TNF Receptor II(green) in HeLa using TNF Receptor II antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using TNF Receptor II antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.