

CD147 Rabbit mAb


ABZ5231

Product Overview


Name	CD147 Rabbit mAb
Accession(Primary)	P18572
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse
Calculated MW	Calculated MW: 42 kDa; Observed MW: 38-58 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:1000-1:5000; ICC/IF-1:100-1:200; FC-1:100-1:200; IP-1:20-1:50

Target information (P18572)


Gene Name	Bsg
About Gene Symbol	



Immunohistochemistry analysis of paraffin-embedded mouse lung using CD147 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of CD147 in mouse liver lysates using CD147 antibody.



Western blot analysis of CD147/ Basigin in HeLa,
C2C12 lysates using CD147/ Basigin antibody